

## 5.46 (NMR Spectroscopy and Organic Structure Determination) Problem Sets

Steven Labalme

February 23, 2025

# Contents


<b>1</b>	<b>Working with NMR Spectra</b>	<b>1</b>
<b>2</b>	<b>Signal-To-Noise Ratio</b>	<b>4</b>

## 1 Working with NMR Spectra


2/19: 1. One of the directories in the 5.46\_NMR\_data folder is called 'glucose2\_320mM'; inside this directory are a number of datasets which were acquired at different times or using different NMR experiments.

- a) This dataset was, surprisingly, generated from a 320 mM sample of glucose. Glucose exists as a mixture of  $\alpha$ - and  $\beta$ -anomers. Draw the structures of each in chair form, and label each anomer.

*Answer.*



(a)  $\alpha$ -D-glucose.




(b)  $\beta$ -D-glucose.

□

- b) Find experiment #4, load it into MestReNova or TopSpin, and generate a properly phased spectrum. Use the software to integrate the peaks or regions, and expand the spectrum to include only the region around water to include the anomeric protons.

*Answer.*



□

- c) Measure the  $^3J_{\text{HH}}$  coupling constant and the  $^1J_{\text{CH}}$  coupling constant for each anomeric proton. Make a table showing the chemical shift for each anomeric proton, the relative integrals (which should total 100%), and the two coupling constants.

*Answer.*

Pyranose	Chemical Shift	Integral	$^3J_{\text{HH}}$	$^1J_{\text{CH}}$
$\alpha$ -D-glucose	5.24 ppm	61%	3.81 Hz	169.50 Hz
$\beta$ -D-glucose	4.65 ppm	39%	7.96 Hz	164.62 Hz

□

- d) Using the fact that the  $\beta$ -anomer is expected to be about 60% of the total, rationalize the variation in chemical shift and coupling constants for each.

*Answer.* Contrary to the expected integrals, the peak characteristic of the  $\alpha$ -anomer integrates to approximately 61% of the total in the provided data file. Perhaps the anomers were taken out of equilibrium before the spectrum was taken and are in the process of reequilibrating.

Regardless, the  $\beta$ -anomer's anomeric proton should be more shielded because of the anomeric effect's no-bond resonance placing a  $\delta^-$  on that proton; such resonance does not occur in the  $\alpha$ -anomer. This effect is clearly observed, as the  $\beta$ -anomer's doublet is 0.59 ppm upfield from the  $\alpha$ -anomer's doublet. Note that the downfield doublets were originally identified as the anomeric peaks because the anomeric protons should be the most deshielded due to their position on the glucose hemiacetal (the position nearest to the most electron-withdrawing hydroxyl groups).

Additionally, the  $\beta$ -anomer's  $^3J_{\text{HH}}$  coupling constant should be larger than the  $\alpha$ -anomer's due to the Karplus correlation. Indeed, the  $\beta$ -anomer's anomeric proton is oriented antiperiplanar to the C2 proton, while the  $\alpha$ -anomer's anomeric proton is oriented gauche. This effect is also clearly observed, as the  $\beta$ -anomer has a  $^3J_{\text{HH}}$  value 4.15 Hz greater than the  $\alpha$ -anomer's  $^3J_{\text{HH}}$  value.

Lastly, the  $\beta$ -anomer's  $^1J_{\text{CH}}$  coupling constant should be smaller than the  $\alpha$ -anomer's because the relevant C–H bond is longer in the  $\beta$ -anomer. This bond elongation is again a result of the anomeric effect, specifically  $n_{\text{O}} \rightarrow \sigma_{\text{CH}}^*$  donation. This effect is again observed, as the  $\beta$  anomer has a  $^1J_{\text{CH}}$  value 4.88 Hz smaller than the  $\alpha$ -anomer's  $^1J_{\text{CH}}$  value. □

2. A second directory in the data folder is called 'glucose1.320mM'; this sample is also glucose at 320 mM with one difference – this glucose is fully  $^{13}\text{C}$  labeled. Analyze the anomeric region of spectrum #4 in the same way as you did for the first sample. Explain the difference in apparent anomeric ratio as best you can.

*Answer.*

Pyranose	Chemical Shift	Integral	$^3J_{\text{HH}}$	$^1J_{\text{CH}}$
$\alpha$ -D-glucose	5.24 ppm	62%	??	169.92 Hz
$\beta$ -D-glucose	4.65 ppm	39%	??	165.08 Hz

Of note is that the central chemical shift values do not change (as expected), and the coupling constants only increase marginally (perhaps insignificantly). The  $^1J_{\text{CH}}$  values are also difficult to measure — especially for the  $\alpha$ -anomer — because of the appearance of apparently new couplings that split the peaks into pseudo-quartets.

Additionally, there does not appear to be a significant difference in the anomeric ratio. Rather, the  $\alpha : \beta$  ratio appears to be approximately 62 : 38 instead of 61 : 39. □

3. For the spectra you have been working with, measure the amount of water relative to the glucose. Estimate the amount of water, and think about where this water signal comes from. Which glucose sample has more water in it?

*Answer.* The residual solvent peak at 4.79 ppm in glucose2\_320mM.4.fid integrates to 3.82 times the height of the sum of the anomeric peaks. This means that for each glucose anomeric proton in solution, there are 3.82 water protons. Since there is one glucose anomeric proton per glucose molecule but two water protons per molecule, it follows that there are 1.91 water molecules for each glucose molecule in solution. Thus, the concentration of water in this sample is

$$1.91 \times 320 \text{ mM} = 611 \text{ mM}$$

By an analogous calculation, the concentration of water in glucose1\_320mM.4.fid is

$$2.01 \times 320 \text{ mM} = 643 \text{ mM}$$

Thus, glucose1\_320mM.4.fid has more water in it. The water signal most likely comes partially from leftover  $\text{H}_2\text{O}$  contamination in the isotopically enriched  $\text{D}_2\text{O}$  solvent, and partially from proton-deuterium exchange between  $\text{D}_2\text{O}$  and undeuterated glucose.  $\square$

## 2 Signal-To-Noise Ratio

2/25: We often measure the sensitivity of a system using a standard sample and a standard experimental design; this allows us to compare one NMR system to another, or one probe to another in the same NMR system. It also allows us to estimate the result we can expect when we run a real sample (one that we care about).

For this exercise we will use data collected for adenosine in d6-DMSO to compare some of the systems in the DCIF with respect to both proton and carbon acquisition. The method is as follows — choose the largest adenosine signal in the spectrum (do not use the water or DMSO peaks) and use either MestReNova (the non-interactive approach is fine) or Topspin to measure SNR (the signal-to-noise ratio) using the software. For the proton spectrum one of the aromatic singlets will be largest; for the carbon spectrum, it will be one of the sugar carbons.

1. For the 20 mM adenosine data acquired at 500 MHz, experiments 1 and 2 are proton acquisitions using 16 and 128 scans, respectively. Measure the SNR for each and confirm that the increase in SNR with scans matches what you expect.
2. Compare experiment 1 of the 20 mM adenosine sample from the 500 with experiment 1 of the 10 mM sample from the 500 (labeled 5.46\_2025\_adenosine\_10mM\_500), and again confirm that the difference in SNR is what you expect.
3. Now compare the SNR for the 10 mM adenosine samples for both proton (experiment 1) and carbon (experiment 2) across all three NMR spectrometers (400, 500, and 600). For both proton and carbon, how many scans would you need at 400 MHz and 500 MHz to give you the same SNR as you see at 600 MHz?
4. From the results in Q3, use the [DCIF fee schedule](#) to calculate the cost of achieving the same SNR for the AVIII 400, Neo500, and Neo600. Assume the proton spectra take 6 min and the carbon spectra take 72 min.
5. Finally, using the 10 mM 400 MHz proton and carbon spectra, measure the SNR with different line broadening (apodization) values — 0 Hz, 1 Hz, 5 Hz, 10 Hz, and 20 Hz. Which line broadening gives the highest SNR for each nucleus? Why?